

## MemoPart™ VSD Occluder Instruction for Use



### **Shanghai Shape Memory Alloy Co., Ltd.**

Address: 1F and 5F, Tower 41, No. 258 XinZhuan Road, Songjiang High-Tech Park,  
CaoHeJing Development District, 201612 Shanghai, P.R. China

Code: 201612

Tel : +86-21-37013390

Fax : +86-21-37013391

Web: <https://en.scientechmed.com/>



### **Lepu Medical (Europe) Cooperatief U.A.**

Address: Abe Lenstra Boulevard 36, 8448 JB, Heerenveen, the Netherlands

Tel: +31-515-573399

Fax: +31-515-760020



**SHSMA**

**Shanghai Shape Memory Alloy Co., Ltd.**

## TABLE OF CONTENT

1.	Product Description.....	1
2.	Indications for Use .....	1
3.	Contraindications .....	1
4.	Intended Patient Populations .....	2
5.	Intended Users.....	2
6.	Warning.....	2
7.	Precautions .....	2
8.	Adverse Events.....	3
9.	Clinical Benefits to be expected .....	3
10.	Summary of Safety and Clinical Performance (SSCP) .....	3
11.	Materials recommended for use with this device .....	3
12.	Directions for Use .....	3
13.	Disposal.....	4
14.	Storage and Expiry Date .....	4
15.	Symbol Definitions .....	4
16.	DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY .....	4
17.	Specification and Recommend sheath size .....	4

**Read “Instruction for Use” carefully before use, Users should have the interventional treating practice or be guided by professionals.**

- Use on or before the last day of the expiration month noted on the product packaging.
- This device is sterilized using ethylene oxide and is for single use only. Do not reuse or resterilize. Attempts to resterilize the device may result in device malfunction, inadequate sterilization, or patient harm.
- Do not use the device if the packaging sterile barrier is open or damaged.

### **1. Product Description**

MemoPart™ VSD Occluder is a self-expanding double-disc nitinol mesh occlusion device. The 2 discs are connected by a short waist which corresponds to the defect size. Polyester fabric is securely sewn to each disc to secure the occlusion. The device is visible under X-ray.

### **2. Intended purpose**

MemoPart™ VSD Occluder is specially designed for transcatheter closure of congenital ventricular septal defect.

#### **2.1 Indications for Use**

- *Membranous VSD Occluder*  
The MemoPart™ Membranous VSD Occluder is used for minimally invasive transcatheter closure of perimembranous ventricular septal defects.
- *Muscular VSD Occluder*
  - 1) Patient with a complex ventricular septal defect (VSD) of significant size to warrant closure (large volume left-to right shunt, pulmonary hypertension, and/or clinical symptoms of congestive heart failure) ;
  - 2) Patients who are considered to be at high risk for standard transatrial or transarterial surgical closure based on anatomical conditions and/or based on overall medical condition. High-risk anatomical factors for transatrial or transarterial surgical closure include patients:
    - Requiring left ventriculotomy or an extensive right ventriculotomy;
    - With a failed previous VSD closure;
    - With multiple apical and/or anterior muscular VSDs (“Swiss cheese septum”);
    - With posterior apical VSDs covered by trabeculae.

#### **2.2 Contraindications**

- *For membranous ventricular septal defect:*
  - ♦ Presence of thrombi in the chambers of the heart and hemorrhagic diseases such as active ulcer;
  - ♦ Recent infection is found;

- ♦ Patients with very small vessels which are inadequate to accommodate the appropriate sheath size;
  - ♦ The patient with a distance from the edge of the defect to the aorta or to the tricuspid valve less than 3mm;
  - ♦ Anatomy in which the Membranous VSD Occluder would interfere with the aortic or atrioventricular valves;
  - ♦ Patients with coagulation disorder who are unable to take antiplatelet or anticoagulant therapy;
  - ♦ Patients with an intra-cardiac mass or vegetation.
- *For muscular ventricular septal defect:*
- ♦ Patients who do not have a Muscular VSD or if the Muscular VSD is from a heart attack;
  - ♦ Patients with an infection anywhere in the body recently. They may receive the device only after the infection is gone;
  - ♦ Patients unable to take aspirin (unless they can take other anti-platelet agents for 6 months);
  - ♦ Patient whose heart or veins are very small, or if he cannot take undergo the procedure;
  - ♦ Patients with defects less than 4 mm distance from the semilunar (aortic and pulmonary) and atrioventricular valves (mitral and tricuspid);
  - ♦ Patients with severely increased pulmonary vascular resistance above 7 Wood units and a right-to-left shunt and documented irreversible pulmonary vascular disease.

### 3. Intended Patient Populations

- 1) This device is suitable for any patients with the above disease or medical condition except for whom with contraindications. The safety and efficacy of MemoPart™ VSD Occluder in pregnant women or men intending to father children, nursing mothers and the immunocompromised patients have not been established.
- 2) As for weight of intended patients:
  - MemoPart™ VSD Occluder is recommended in the transcatheter closure of VSD in patients weighting  $\geq 10\text{kg}$ ;
  - For patients weighting  $< 10\text{kg}$ , MemoPart™ VSD Occluder should be used under the suggestion of clinical professional.

### 4. Intended Users

The MemoPart™ VSD Occluder should only be used by a physician who are trained in transcatheter defect closure techniques. The physician should determine which patients are candidates for procedures that use this device.

Before use, the operator should have a full understanding of the using instruction, warnings and precautions.

### 5. Warning

- 6.1 Patients who are allergic to nickel may have an allergic reaction to this device.
- 6.2 This device should only be used by physicians who have been trained in transcatheter techniques and who should determine which patients are suitable candidates for procedures using this device.
- 6.3 Physicians must be prepared to deal with urgent situations, such as device embolization, which require removal of the device. This includes the availability of an on-site surgeon.
- 6.4 This product is fit for magnetic resonance imaging.
- 6.5 There is a variety of pocket-type (aneurysm of membranous septum) ventricular septal defect and incomplete occlusion may take septal defects with big entrances and multiple exits, so appropriate occluder should be chosen with care. Angiography in left ventricle should be performed before releasing the occluder to determine if the shunt is present or not.
- 6.6 It is important to choose the appropriate sheath to the ventricular septal defect occluder (see Table 1); if the sheath cannot match the occluder, it would be difficult to push or retrieve the occluder which may even cause damages to the blood vessels.

### 6. Precautions

Before foreign like product happened the unexpected situations during the using were the cardiac tamponade caused by cardiac wall perforation, the dropping of occluder, auriculo-ventricular block, aortic insufficiency, tricuspid incompetence and hemolysis. Therefore, it should be pay an attention to prevention in the operation. And it's required that the institute which has the conditions for cardiothoracic surgery could use this occluder. In case the cardiac tamponade caused by cardiac wall perforation was occurred, the pericardiocentesis drainage should be immediately implemented, if the occluder was dropped after released, the emergency surgery should be implemented to remove the occluder. MR Conditional: Non-clinical testing has demonstrated that the VSD occluder is MR Conditional.

Patients can be scanned safely immediately after implantation under the following conditions:

- Static magnetic field of 1.5T and 3.0T;
- Maximum allowable spatial gradient of the magnetic field is 20T/m in 3.0T and 40T/m in 1.5T MR system;
- Maximum whole-body specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning in normal operating mode.

- The presence of this implant might produce an image artifact.

## 7. Adverse Events

Potential adverse events may occur during or after a procedure placing this device may include, but are not limited to:

Death, Residual shunt, Cardiac arrhythmia, Heart block/Conductive block, Thrombosis, Dislocation/malposition/device migration, Device embolization, Infection including endocarditis, Fever, Cardiac perforation/tamponade, Cardiac erosion, Pericardial effusion, Pleural effusion, Valvular insufficiency/regurgitation, Valve damage, Vessel trauma/damage, Hematoma, Bleeding/Blood transfusion, Hemolysis, Arterio-venous fistula, Hypertension/hypotension, CVA/Stroke/TIA, Subaortic Stenosis, Ventricular outflow tract stenosis, Chest pain, Anesthesia reactions, Headache/migraine, Pulmonary hypertension, Cyanosis, Vomiting, Pulmonary edema, Arterial pulse loss, Brachial plexus injury, Device fracture, Syncope, Atelectasis, Cardiomyopathy, Abdominal gassiness, Stridor, Allergic reaction.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 8. Clinical Benefits to be expected

- High procedural success;
- High complete closure rate;
- Low rate of clinical adverse event.

## 9. Summary of Safety and Clinical Performance (SSCP)

The summary of safety and performance for the device is available on the website of EUDAMED: <https://ec.europa.eu/tools/eudamed>.

## 10. Materials recommended for use with this device

- ♦ 0.035-inch guidewire of sufficient stiffness
- ♦ 20 ml luer-lock syringe

## 11. Directions for Use

- 12.1 Perform preoperative disinfection and spread out the operating towel and drape.
- 12.2 Local anesthesia or general anesthesia.
- 12.3 Seldinger puncture at right femoral vein and artery is performed.
- 12.4 Place the leakage-proof sheath.
- 12.5 Inject 100μ/kg of heparin into the vein which would be increased by 1000-2000μ every hour.

- 12.6 Perform routine cardiac catheter examination and determine pressure inside each cardiac chamber as well as the blood oxygen content.
- 12.7 Deliver the angiography catheter of right coronary into the left ventricle and then into right ventricle via the ventricular septal defect; deliver the guiding wire through the catheter into precava or pulmonary artery (see Fig 2).
- 12.8 Deliver the right cardiac catheter into the pulmonary artery via the vein and into the snare through the catheter; hitch the guiding wire and pull it out of the body after which the artery-ventricular septal defect-vein track has been established.




Fig. 2 Operation chart for artery-vein track

- 12.9 Insert the guiding wire into the delivery sheath and deliver it into the left ventricle. Choose the occluder with appropriate size according to the results of angiography; mostly the disc diameter of occluder would be 2mm larger than the determined diameter.
- 12.10 Put the occluder into heparinized saline, connect it with the pusher, rotate it clockwise for 4 rounds and pull it into the short sheath. Push in heparinized saline rapidly through the connecting tube of the haemostatic by-pass valve and drain all the bubbles out of the occlusion device.
- 12.11 Insert the dilator into the longer sheath and push forward under X-ray; open the disc at the side of the left ventricle and retrieve it back to ventricular septum; if the encountered resistance is tested to be around the ventricular septum by echocardiography, then fix the pusher and retrieve the sheath as well as release the disc of the occluder at the side of the right ventricle. Push and pull the occluder to see if it is fixed; if echocardiography shows that it would not impact the function of aortic valve and tricuspid valve and the angiography shows that it would not impact the function of aortic valve, and then the pushing pole could be rotated counterclockwise to release the occluder.




Fig. 3 Operation chart for VSD

- 12.12 During the operation, TTE is performed from the five-chamber view at the

cardiac apex and the short-axis view. The occluder could only be released when all results from these two views have shown that the occluder is closely next to the residual edge of the ventricular septum.

12.13 Push out the sheath and finish the operation.

12.14 Take orally 3-5 mg/kg aspirin after operation for 6 months.

## 12. Disposal

Treat all disposable devices appropriately according to the local requirements for medical device waste disposal.

## 13. Storage and Expiry Date

14.1 Store in a dry place without ultraviolet radiation and should be stored in a well-ventilated room with a relative humidity of no more than 80% and no corrosion gas.

14.2 The period of validity is five years when this product is stored under specified condition.

## 14. Symbol Definitions

SYMBOL	EXPLANATION OF SYMBOL
	MANUFACTURER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	DATE OF MANUFACTURE
	USE-BY DATE
	BATCH CODE
	CATALOGUE NUMBER
	SERIAL NUMBER
	STERILIZED USING ETHYLENE OXIDE
	DOUBLE STERILE BARRIER SYSTEM
	SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE
	SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE AND STERILIZED USING ETHYLENE OXIDE

	DO NOT RESTERILIZE
	DO NOT RE-USE
	DO NOT USE IF PACKAGE IS DAMAGED
	CONSULT INSTRUCTIONS FOR USE
	MEDICAL DEVICE
	UNIQUE DEVICE IDENTIFIER
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	MR CONDITIONAL
	CE MARKING AND IDENTIFICATION NUMBER OF NOTIFIED BODY
The color of EO Sterilization indicator on the package bag turns yellow after EO sterilization.	

## 15. DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

Descriptions or specifications in Shanghai Shape Memory Alloy Co., Ltd. printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Shanghai Shape Memory Alloy Co., Ltd will not be responsible for any direct, incidental, or consequential damages resulting from the abnormal use of the product.

## 16. Specification and Recommend sheath size

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
SQFDQ- I a04	04	8.0±1.0	5.0±1.5	4.0±1.0	8.0±1.0	7-8F	Muscular defect
SQFDQ- I a05	05	9.0±1.0	5.0±1.5	5.0±1.0	9.0±1.0	7-8F	
SQFDQ- I a06	06	10.0±1.0	5.0±1.5	6.0±1.0	10.0±1.0	8-9F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
SQFDQ- I a07	07	11.0±1.0	5.0±1.5	7.0±1.0	11.0±1.0	8-9F	
SQFDQ- I a08	08	12.0±1.0	5.0±1.5	8.0±1.0	12.0±1.0	8-9F	
SQFDQ- I a09	09	13.0±1.0	5.0±1.8	9.0±1.0	13.0±1.0	8-9F	
SQFDQ- I a10	10	14.0±1.0	5.0±1.8	10.0±1.2	14.0±1.0	9-10F	
SQFDQ- I a12	12	16.0±1.0	5.0±1.8	12.0±1.2	16.0±1.0	9-10F	
SQFDQ- I a14	14	18.0±1.0	5.0±1.8	14.0±1.5	18.0±1.0	10-12F	
SQFDQ- I a16	16	20.0±1.0	5.0±1.8	16.0±1.5	20.0±1.0	10-12F	
SQFDQ- I a18	18	22.0±1.0	5.0±1.8	18.0±1.5	22.0±1.0	10-12F	
SQFDQ- I b04	04	10.0±1.0	7.0±1.5	4.0±1.0	8.0±1.0	7-8F	
SQFDQ- I b05	05	11.0±1.0	7.0±1.5	5.0±1.0	9.0±1.0	7-8F	
SQFDQ- I b06	06	12.0±1.0	7.0±1.5	6.0±1.0	10.0±1.0	8-9F	
SQFDQ- I b07	07	13.0±1.0	7.0±1.5	7.0±1.0	11.0±1.0	8-9F	
SQFDQ- I b08	08	14.0±1.0	7.0±1.5	8.0±1.0	12.0±1.0	8-9F	Muscular defect
SQFDQ- I b09	09	15.0±1.0	7.0±1.8	9.0±1.0	13.0±1.0	8-9F	
SQFDQ- I b10	10	16.0±1.0	7.0±1.8	10.0±1.2	14.0±1.0	9-10F	
SQFDQ- I b12	12	18.0±1.0	7.0±1.8	12.0±1.2	16.0±1.0	9-10F	
SQFDQ- I b14	14	20.0±1.0	7.0±1.8	14.0±1.5	18.0±1.0	10-12F	
SQFDQ- I b16	16	22.0±1.0	7.0±1.8	16.0±1.5	20.0±1.0	10-12F	
SQFDQ- I b18	18	24.0±1.0	7.0±1.8	18.0±1.5	22.0±1.0	10-12F	
SQFDQ- I c04	04	14.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	
SQFDQ- I c05	05	15.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
SQFDQ- I c06	06	16.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
SQFDQ- I c07	07	17.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
SQFDQ- I c08	08	18.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
SQFDQ- I c09	09	19.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	
SQFDQ- I c10	10	20.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
SQFDQ- I c12	12	22.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
SQFDQ- I c14	14	24.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	Muscular defect
SQFDQ- I c16	16	26.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	
SQFDQ- I c18	18	28.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	
SQFDQ- I d04	04	18.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	
SQFDQ- I d05	05	19.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
SQFDQ- I d06	06	20.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
SQFDQ- I d07	07	21.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
SQFDQ- I d08	08	22.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
SQFDQ- I d09	09	23.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	
SQFDQ- I d10	10	24.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
SQFDQ- I d12	12	26.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
SQFDQ- I d14	14	28.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	
SQFDQ- I d16	16	30.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	
SQFDQ- I d18	18	32.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	
WTSQFDQ- I a04	04	8.0±1.0	5.0±1.5	4.0±1.0	8.0±1.0	7-8F	Muscular defect

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
WTSQFDQ- I a05	05	9.0±1.0	5.0±1.5	5.0±1.0	9.0±1.0	7-8F	
WTSQFDQ- I a06	06	10.0±1.0	5.0±1.5	6.0±1.0	10.0±1.0	8-9F	
WTSQFDQ- I a07	07	11.0±1.0	5.0±1.5	7.0±1.0	11.0±1.0	8-9F	
WTSQFDQ- I a08	08	12.0±1.0	5.0±1.5	8.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I a09	09	13.0±1.0	5.0±1.8	9.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I a10	10	14.0±1.0	5.0±1.8	10.0±1.2	14.0±1.0	9-10F	
WTSQFDQ- I a12	12	16.0±1.0	5.0±1.8	12.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I a14	14	18.0±1.0	5.0±1.8	14.0±1.5	18.0±1.0	10-12F	
WTSQFDQ- I a16	16	20.0±1.0	5.0±1.8	16.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I a18	18	22.0±1.0	5.0±1.8	18.0±1.5	22.0±1.0	10-12F	
WTSQFDQ- I b04	04	10.0±1.0	7.0±1.5	4.0±1.0	8.0±1.0	7-8F	
WTSQFDQ- I b05	05	11.0±1.0	7.0±1.5	5.0±1.0	9.0±1.0	7-8F	
WTSQFDQ- I b06	06	12.0±1.0	7.0±1.5	6.0±1.0	10.0±1.0	8-9F	Muscular defect
WTSQFDQ- I b07	07	13.0±1.0	7.0±1.5	7.0±1.0	11.0±1.0	8-9F	
WTSQFDQ- I b08	08	14.0±1.0	7.0±1.5	8.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I b09	09	15.0±1.0	7.0±1.8	9.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I b10	10	16.0±1.0	7.0±1.8	10.0±1.2	14.0±1.0	9-10F	
WTSQFDQ- I b12	12	18.0±1.0	7.0±1.8	12.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I b14	14	20.0±1.0	7.0±1.8	14.0±1.5	18.0±1.0	10-12F	
WTSQFDQ- I b16	16	22.0±1.0	7.0±1.8	16.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I b18	18	24.0±1.0	7.0±1.8	18.0±1.5	22.0±1.0	10-12F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
WTSQFDQ- I c04	04	14.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	
WTSQFDQ- I c05	05	15.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
WTSQFDQ- I c06	06	16.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I c07	07	17.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I c08	08	18.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
WTSQFDQ- I c09	09	19.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	
WTSQFDQ- I c10	10	20.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I c12	12	22.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
WTSQFDQ- I c14	14	24.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I c16	16	26.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	
WTSQFDQ- I c18	18	28.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	
WTSQFDQ- I d04	04	18.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	Muscular defect
WTSQFDQ- I d05	05	19.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
WTSQFDQ- I d06	06	20.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I d07	07	21.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I d08	08	22.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
WTSQFDQ- I d09	09	23.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	
WTSQFDQ- I d10	10	24.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I d12	12	26.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
WTSQFDQ- I d14	14	28.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I d16	16	30.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
WTSQFDQ- I d18	18	32.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	Muscular defect
SQFDQ-IIa04	04	8.0±0.8	1.8±0.5	4.0±0.8	8.0±0.8	6-7F	Membranous defect
SQFDQ-IIa05	05	9.0±0.8	1.8±0.5	5.0±0.8	9.0±0.8	6-7F	
SQFDQ-IIa06	06	10.0±0.8	1.8±0.5	6.0±0.8	10.0±0.8	7-8F	
SQFDQ-IIa07	07	11.0±0.8	1.8±0.5	7.0±0.8	11.0±0.8	7-8F	
SQFDQ-IIa08	08	12.0±0.8	1.8±0.5	8.0±0.8	12.0±0.8	7-8F	
SQFDQ-IIa09	09	13.0±0.8	1.8±0.5	9.0±0.8	13.0±0.8	8-9F	
SQFDQ-IIa10	10	14.0±0.8	1.8±0.5	10.0±0.8	14.0±0.8	8-9F	
SQFDQ-IIa12	12	16.0±0.8	1.8±0.5	12.0±0.8	16.0±0.8	9-10F	
SQFDQ-IIa14	14	18.0±0.8	1.8±0.5	14.0±0.8	18.0±0.8	9-10F	
SQFDQ-IIa16	16	20.0±0.8	1.8±0.5	16.0±0.8	20.0±0.8	10-12F	
SQFDQ-IIa18	18	24.0±0.8	1.8±0.5	18.0±0.8	22.0±0.8	10-12F	
SQFDQ-IIa20	20	26.0±0.8	1.8±0.5	20.0±0.8	24.0±0.8	12-14F	
SQFDQ-IIb04	04	8.0±0.8	3.5±1.0	4.0±0.8	8.0±0.8	6-7F	Membranous defect
SQFDQ-IIb05	05	9.0±0.8	4.0±1.0	5.0±0.8	9.0±0.8	6-7F	
SQFDQ-IIb06	06	10.0±0.8	4.0±1.0	6.0±1.0	10.0±0.8	7-8F	
SQFDQ-IIb07	07	11.0±1.0	4.0±1.0	7.0±1.0	11.0±1.0	7-8F	
SQFDQ-IIb08	08	12.0±1.0	4.0±1.0	8.0±1.0	12.0±1.0	7-8F	
SQFDQ-IIb09	09	13.0±1.0	4.5±1.0	9.0±1.2	13.0±1.0	8-9F	
SQFDQ-IIb10	10	14.0±1.5	4.5±1.0	10.0±1.2	14.0±1.5	8-9F	
SQFDQ-IIb12	12	16.0±1.5	4.5±1.0	12.0±1.5	15.0±1.5	9-10F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
SQFDQ-IIb14	14	18.0±1.5	4.5±1.0	14.0±1.5	17.0±1.5	9-10F	Membranous defect
SQFDQ-IIb16	16	22.0±1.5	5.0±1.0	16.0±1.5	20.0±1.5	10-12F	
SQFDQ-IIb18	18	24.0±1.5	5.0±1.0	18.0±1.8	22.0±1.5	10-12F	
SQFDQ-IIb20	20	26.0±1.5	5.0±1.0	20.0±1.8	24.0±1.5	12-14F	
WTSQFDQ-IIa04	04	8.0±0.8	1.8±0.5	4.0±0.8	8.0±0.8	6-7F	
WTSQFDQ-IIa05	05	9.0±0.8	1.8±0.5	5.0±0.8	9.0±0.8	6-7F	
WTSQFDQ-IIa06	06	10.0±0.8	1.8±0.5	6.0±0.8	10.0±0.8	7-8F	
WTSQFDQ-IIa07	07	11.0±0.8	1.8±0.5	7.0±0.8	11.0±0.8	7-8F	
WTSQFDQ-IIa08	08	12.0±0.8	1.8±0.5	8.0±0.8	12.0±0.8	7-8F	
WTSQFDQ-IIa09	09	13.0±0.8	1.8±0.5	9.0±0.8	13.0±0.8	8-9F	
WTSQFDQ-IIa10	10	14.0±0.8	1.8±0.5	10.0±0.8	14.0±0.8	8-9F	
WTSQFDQ-IIa12	12	16.0±0.8	1.8±0.5	12.0±0.8	16.0±0.8	9-10F	
WTSQFDQ-IIa14	14	18.0±0.8	1.8±0.5	14.0±0.8	18.0±0.8	9-10F	Membranous defect
WTSQFDQ-IIa16	16	20.0±0.8	1.8±0.5	16.0±0.8	20.0±0.8	10-12F	
WTSQFDQ-IIa18	18	24.0±0.8	1.8±0.5	18.0±0.8	22.0±0.8	10-12F	
WTSQFDQ-IIa20	20	26.0±0.8	1.8±0.5	20.0±0.8	24.0±0.8	12-14F	
WTSQFDQ-IIb04	04	8.0±0.8	3.5±1.0	4.0±0.8	8.0±0.8	6-7F	
WTSQFDQ-IIb05	05	9.0±0.8	4.0±1.0	5.0±0.8	9.0±0.8	6-7F	
WTSQFDQ-IIb06	06	10.0±0.8	4.0±1.0	6.0±1.0	10.0±0.8	7-8F	
WTSQFDQ-IIb07	07	11.0±1.0	4.0±1.0	7.0±1.0	11.0±1.0	7-8F	
WTSQFDQ-IIb08	08	12.0±1.0	4.0±1.0	8.0±1.0	12.0±1.0	7-8F	



Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
WTSQFDQ-IIb09	09	13.0±1.0	4.5±1.0	9.0±1.2	13.0±1.0	8-9F	
WTSQFDQ-IIb10	10	14.0±1.5	4.5±1.0	10.0±1.2	14.0±1.5	8-9F	
WTSQFDQ-IIb12	12	16.0±1.5	4.5±1.0	12.0±1.5	15.0±1.5	9-10F	
WTSQFDQ-IIb14	14	18.0±1.5	4.5±1.0	14.0±1.5	17.0±1.5	9-10F	Membranous defect
WTSQFDQ-IIb16	16	22.0±1.5	5.0±1.0	16.0±1.5	20.0±1.5	10-12F	
WTSQFDQ-IIb18	18	24.0±1.5	5.0±1.0	18.0±1.8	22.0±1.5	10-12F	
WTSQFDQ-IIb20	20	26.0±1.5	5.0±1.0	20.0±1.8	24.0±1.5	12-14F	
SQFDQ-III04	04	12.0±1.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
SQFDQ-III05	05	13.0±1.0	4.0±1.5	5.0±1.0	9.0±1.0	8-9F	
SQFDQ-III06	06	14.0±1.0	4.0±1.5	6.0±1.0	10.0±1.0	8-9F	
SQFDQ-III07	07	15.0±1.0	4.0±1.5	7.0±1.2	11.0±1.0	8-9F	
SQFDQ-III08	08	16.0±1.2	4.0±1.5	8.0±1.2	12.0±1.2	9-10F	
SQFDQ-III09	09	17.0±1.2	4.5±1.5	9.0±1.2	13.0±1.2	9-10F	
SQFDQ-III10	10	18.0±1.2	4.5±1.5	10.0±1.5	14.0±1.2	9-10F	
SQFDQ-III12	12	20.0±1.5	4.5±1.5	12.0±1.5	16.0±1.2	10-12F	
SQFDQ-III14	14	22.0±1.5	4.5±1.5	14.0±1.8	18.0±1.5	10-12F	
SQFDQ-III16	16	24.0±1.5	5.0±1.5	16.0±1.8	20.0±1.5	10-12F	
SQFDQ-III18	18	26.0±1.5	5.0±1.5	18.0±1.8	22.0±1.5	12-14F	
WTSQFDQ-III04	04	12.0±1.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	Membranous defect
WTSQFDQ-III05	05	13.0±1.0	4.0±1.5	5.0±1.0	9.0±1.0	8-9F	
WTSQFDQ-III06	06	14.0±1.0	4.0±1.5	6.0±1.0	10.0±1.0	8-9F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
WTSQFDQ-III07	07	15.0±1.0	4.0±1.5	7.0±1.2	11.0±1.0	8-9F	
WTSQFDQ-III08	08	16.0±1.2	4.0±1.5	8.0±1.2	12.0±1.2	9-10F	
WTSQFDQ-III09	09	17.0±1.2	4.5±1.5	9.0±1.2	13.0±1.2	9-10F	
WTSQFDQ-III10	10	18.0±1.2	4.5±1.5	10.0±1.5	14.0±1.2	9-10F	
WTSQFDQ-III12	12	20.0±1.5	4.5±1.5	12.0±1.5	16.0±1.2	10-12F	
WTSQFDQ-III14	14	22.0±1.5	4.5±1.5	14.0±1.8	18.0±1.5	10-12F	
WTSQFDQ-III16	16	24.0±1.5	5.0±1.5	16.0±1.8	20.0±1.5	10-12F	
WTSQFDQ-III18	18	26.0±1.5	5.0±1.5	18.0±1.8	22.0±1.5	12-14F	
SQFDQ-IV04	04	9.0±2.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
SQFDQ-IV05	05	10.0±2.0	3.5±1.5	5.0±0.8	9.0±1.0	7-8F	
SQFDQ-IV06	06	11.0±2.0	4.0±1.5	6.0±1.0	10.0±1.0	7-8F	
SQFDQ-IV07	07	12.0±2.5	4.0±1.5	7.0±1.0	11.0±1.0	8-9F	
SQFDQ-IV08	08	13.0±2.5	4.5±1.5	8.0±1.2	12.0±1.2	8-9F	Membranous defect
SQFDQ-IV09	09	14.0±2.5	5.0±1.5	9.0±1.2	13.0±1.2	9-10F	
SQFDQ-IV10	10	17.0±2.5	5.0±1.5	10.0±1.7	15.0±1.2	9-10F	
SQFDQ-IV12	12	20.0±3.5	5.0±1.5	12.0±1.7	18.0±1.2	10-12F	
SQFDQ-IV14	14	22.0±3.5	5.0±1.5	14.0±1.8	20.0±1.5	10-12F	
SQFDQ-IV16	16	24.0±3.5	5.0±1.5	16.0±1.8	22.0±1.5	10-12F	
WTSQFDQ-IV04	04	9.0±2.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
WTSQFDQ-IV05	05	10.0±2.0	3.5±1.5	5.0±0.8	9.0±1.0	7-8F	
WTSQFDQ-IV06	06	11.0±2.0	4.0±1.5	6.0±1.0	10.0±1.0	7-8F	

Catalogue No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Sheath Size	Applicable defect
WTSQFDQ-IV07	07	12.0±2.5	4.0±1.5	7.0±1.0	11.0±1.0	8-9F	
WTSQFDQ-IV08	08	13.0±2.5	4.5±1.5	8.0±1.2	12.0±1.2	8-9F	
WTSQFDQ-IV09	09	14.0±2.5	5.0±1.5	9.0±1.2	13.0±1.2	9-10F	
WTSQFDQ-IV10	10	17.0±2.5	5.0±1.5	10.0±1.7	15.0±1.2	9-10F	
WTSQFDQ-IV12	12	20.0±3.5	5.0±1.5	12.0±1.7	18.0±1.2	10-12F	
WTSQFDQ-IV14	14	22.0±3.5	5.0±1.5	14.0±1.8	20.0±1.5	10-12F	
WTSQFDQ-IV16	16	24.0±3.5	5.0±1.5	16.0±1.8	22.0±1.5	10-12F	Membranous defect

*Note: The device that is suitable for muscular ventricular septal defect cannot be applied to membranous ventricular septal defect and vice versa.*