

CE 2797

# MemoPart™ VSD Occluder Instruction for Use



## Shanghai Shape Memory Alloy Co., Ltd.

Address: 1F and 5F, Tower 41, No. 258 XinZhuan Road, Songjiang High-Tech Park,  
CaoHeJing Development District, Shanghai 201612, China

Code: 201612

Tel : +86-21-37013390

Fax : +86-21-37013391

Web: <https://en.scientechmed.com/>



## Lepu Medical (Europe) Cooperatief U.A.

Address: Abe Lenstra Boulevard 36, 8448 JB, Heerenveen, the Netherlands

Tel: +31-515-573399

Fax: +31-515-760020



SHSMA

Shanghai Shape Memory Alloy Co., Ltd.

## TABLE OF CONTENT

1.Product Description.....	1
2.Performance characteristics .....	1
3.Intended purpose .....	2
4.Intended Patient Populations .....	3
5.Intended Users .....	3
6.Warning .....	3
7.Precautions .....	4
8.Adverse Events .....	5
9.Expected Clinical Benefits .....	5
10.Summary of Safety and Clinical Performance (SSCP) .....	5
11.Materials recommended for use with this device .....	5
12.Directions for Use .....	5
13.Disposal .....	7
14.Storage and Expiry Date .....	7
15.Symbol Definitions .....	7
16.DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY .....	8
17.Specification and Recommend sheath size .....	9
18.E-IFU.....	17

**Read the “Instruction for Use” carefully before use. Users should possess interventional treatment experience or be guided by professionals.**

- Use on or before the last day of the expiration month noted on the product packaging.
- This device has been sterilized using ethylene oxide and is intended for single use only. Do not reuse or resterilize. Attempts to resterilize the device may result in device malfunction, insufficient sterilization, or patient harm.
- Do not use the device if the sterile packaging barrier is open or damaged.
- An implant card and an implant card information leaflet are supplied with each device. Please fill the implant card after reading the leaflet and provide the card to the patient.

### **1.Product Description**

MemoPart™ VSD Occluder is a self-expanding, double-disc device made from a nickel-titanium alloy wire mesh. The two discs are connected by a short waist corresponding to the size of the ventricular septal defect (VSD). Polyester fabric is securely sewn onto the device to increase the defect closure ability. The device is visible under X-ray.

The company has developed multiple types of occluders to cater to different types of ventricular septal defects, their size, and location, and to reach the desired safety and efficacy levels. The following occluder variants are available:

- ◆ LV disc diameter: 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24,26, 28, 30, and 32 mm.
- ◆ RV disc diameter: 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 20, 22, and 24 mm.
- ◆ Height of connecting waist: 1.8, 3.5, 4, 4.5, 5, 6, 6.5, 7, 7.5, 8, 8.5, 9.5, 10, 10.5, and 12 mm.
- ◆ Connecting waist diameter: 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, and 20 mm.

The corresponding relationship with the catalog number is listed in the table in section 16.

The primary materials used in the device are nitinol (23.8%~68.1%), stainless steel (23.8%~64.2%), and polyester (8.1%~12.0%); all materials used were subjected to biological safety assessment and do not contain CMR and/or endocrine-disrupting substances.

The Occluder Delivery System facilitates the attachment, loading, delivery and deployment of the VSD Occluder.

### **2.Performance characteristics**

This device has the following performance characteristics ensuring its intended clinical performance:

- (1) Nickel-titanium alloy wire constructed the device mesh frame which to facilitate the device deployment;
- (2) Polyester fabric film filled in the device is engineered for the rapid cessation of left-to-right shunts in the ventricular septum upon implantation to prevent ventricular septal defect complications.

(3) Having multiple sizes of LV and RV disc diameter, height of connecting waist and connecting waist diameter, i.e. multiple size of the devices are available for the different defects;

(4) Visible under X-ray is conducive to operation.

### **3.Intended purpose**

MemoPart™ VSD Occluder is specially designed for transcatheter closure of congenital ventricular septal defect.

#### **3.1 Indications for Use**

- For membranous ventricular septal defect:

The MemoPart™ Membranous VSD Occluder is used for minimally invasive transcatheter closure of perimembranous ventricular septal defects.

- For muscular ventricular septal defect:

The Memopart™ Muscular VSD occluder is indicated for the following medical indications:

(1) Patient with a complex ventricular septal defect (VSD) of a significant size that warrants closure (large volume left-to-right shunt, pulmonary hypertension, and/or clinical symptoms of congestive heart failure).

(2) Patients considered to be at high risk for standard transatrial or transarterial surgical closure based on anatomical conditions and/or based on overall medical condition.

High-risk anatomical factors for transatrial or transarterial surgical closure include patients:

- Requiring left ventriculotomy or an extensive right ventriculotomy.
- With a failed previous VSD closure.
- With multiple apical and/or anterior muscular VSDs (“Swiss cheese septum”).
- With posterior apical VSDs covered by trabeculae.

#### **3.2 Contraindications**

- For membranous ventricular septal defect:

- Presence of thrombi in the heart chambers or hemorrhagic diseases such as active ulcer.
- Recent infection.
- Patients with very small vessels that are insufficient to accommodate the appropriate sheath size.
- Patients with a distance from the edge of the defect to the aorta or to the tricuspid valve of less than 3 mm.
- Patients with an anatomy wherein the membranous VSD occluder would interfere with the aortic or atrioventricular valves.
- Patients with coagulation disorder who are unable to take antiplatelet or anticoagulant therapy.

- Patients with an intra-cardiac mass or vegetation.

- For muscular ventricular septal defect:

- Patients who do not have a muscular VSD or if the muscular VSD is from a heart attack.

- Patients with a recent infection anywhere in the body; they may receive the device only after the infection has completely cleared.

- Patients unable to take aspirin (unless they can take other anti-platelet agents for 6 months).

- Patient whose heart or veins are very small, or if they cannot undergo the procedure.

- Patients with defects less than 4 mm distance from the semilunar (aortic and pulmonary) and atrioventricular valves (mitral and tricuspid).

- Patients with severely increased pulmonary vascular resistance above 7 Wood units, a right-to-left shunt, and documented irreversible pulmonary vascular disease.

### **4.Intended Patient Populations**

(1) This device is suitable for any patient with a congenital ventricular septal defect except for those in whom it is contraindicated. The safety and efficacy of MemoPart™ VSD Occluder in pregnant women or men intending to father children, nursing mothers, and the immunocompromised patients have not been established.

(2) Weight-based restrictions of use:

- MemoPart™ VSD Occluder is recommended in the transcatheter closure of VSD in patients weighing  $\geq 10$  kg.

- For patients weighing  $< 10$  kg, MemoPart™ VSD Occluder should only be used after confirmation by a clinical professional.

### **5.Intended Users**

The MemoPart™ VSD Occluder should only be used by a physician who is trained in transcatheter defect closure techniques. The physician should determine the patients who can be candidates for procedures that use this device. Before use, the physician should have a complete understanding of the directions for use, warnings and precautions.

### **6.Warning**

6.1 Patients who are allergic to nickel may have an allergic reaction to this device.

6.2 This device should only be used by physicians who have been trained in transcatheter techniques. The physicians should determine which patients are suitable candidates for procedures using this device.

6.3 Physicians must be prepared to deal with emergency situations, such as device embolization, which require removal of the device. This includes the availability of an on-site surgeon.

6.4 This product is suitable for magnetic resonance imaging.

6.5 The morphology of pocket-type ventricular septal defects, including those that involve a membranous septal aneurysm, exhibit considerable diversity. In situations with cases featuring a large defect entrance and multiple outlets, achieving complete occlusion may not always be feasible pose challenges. Therefore, it is essential to exercise prudence when selecting an appropriate occluder during surgical intervention. Prior to releasing the occluder, a left ventricular angiography should be conducted to confirm the absence of any residual shunt, and only then should the occluder be deployed.

6.6 The device should not be used with delivery systems other than those produced by the manufacturer as the device may be incompatible with the screw thread on the distal end of the pusher of a non-specified delivery system, which may result in technical failures and/or adverse events.

6.7 It is important to choose the appropriate delivery system for the ventricular septal defect occluder (see Table in section 16). If the sheath cannot match the occluder, it would be difficult to push or retrieve the occluder which may cause injury to the blood vessels.

## 7. Precautions

7.1 Unexpected situations during the use of the device includes cardiac tamponade caused by cardiac wall perforation, dropping of the occluder, auriculoventricular block, aortic insufficiency, tricuspid insufficiency, and hemolysis. Therefore, attention should be paid to prevention during the operation. And it is required that the device should be used under institute which has the conditions for cardiothoracic surgery. When cardiac tamponade occurs due to cardiac wall perforation, immediate pericardiocentesis drainage should be carried out. In the event that the occluder is accidentally dropped after being released, emergency surgery must be promptly performed to remove the occluder.

7.2 MR Conditional:

Non-clinical testing has demonstrated that the VSD occluder is MR Conditional. Patients can be safely scanned immediately after implantation under the following conditions:

- Static magnetic field of 1.5T and 3.0T

- Maximum allowable spatial gradient of the magnetic field is 20T/m in 3.0T and 40T/m in 1.5T MR system.
- Maximum whole-body specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning in normal operating mode.
- The presence of this implant may produce an image artifact.

## 8. Adverse Events

Potential adverse events that may occur during or after a procedure for implanting this device may include, but are not limited to:

Death, residual shunt, cardiac arrhythmia, heart block/conductive block, thrombosis, dislocation/malposition/device migration, device embolization, infection including endocarditis, fever, cardiac perforation/tamponade, cardiac erosion, pericardial effusion, pleural effusion, valvular insufficiency/regurgitation, valve damage, vessel trauma/damage, hematoma, bleeding/blood transfusion, hemolysis, arterio-venous fistula, hypertension/hypotension, CVA(Stroke/TIA), subaortic stenosis, ventricular outflow tract stenosis, chest pain, anesthesia reactions, headache/migraine, pulmonary hypertension, cyanosis, vomiting, pulmonary edema, arterial pulse loss, brachial plexus injury, device fracture, syncope, atelectasis, cardiomyopathy, abdominal gassiness, stridor, allergic reaction.

Any serious incident related to the device must be promptly reported to both the manufacturer and the relevant regulatory authority in the member state where the user and/or patient is based.

## 9. Expected Clinical Benefits

- High success rate of device implantation.
- High complete closure rate of ventricular septal defects.
- Low incidence of related clinical adverse events during and after the operation.

## 10. Summary of Safety and Clinical Performance (SSCP)

The summary of safety and performance of the device is available on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>.

Basic UDI-DI: 69330523XZ0002K3

## 11. Materials recommended for use with this device

- 0.035-inch guidewire of sufficient stiffness
- Cardiac catheter
- Delivery system (for matching size, see table in section 16)
- 20 ml luer-lock syringe

## 12. Directions for Use

12.1 Perform preoperative disinfection and spread out the operating towel and

drape.

12.2 Local anesthesia or general anesthesia.

12.3 Perform a Seldinger puncture at right femoral vein and artery.

12.4 Place the leakage-proof sheath.

12.5 Inject 100 µ/kg of heparin into the vein, which should be increased by 1000–2000 µ every hour.

12.6 Perform routine cardiac catheter examination and determine pressure within each cardiac chamber as well as the blood oxygen content.

12.7 Deliver the angiography catheter of right coronary into the left ventricle and then into right ventricle via the ventricular septal defect; deliver the guide wire through the catheter into precava or pulmonary artery (see Fig 2).

12.8 Introduce the right cardiac catheter into the pulmonary artery through the vein and into the snare using the catheter. Capture the guide wire with the snare and carefully withdraw it from the body, establishing the artery-ventricular septal defect-vein track.




Fig. 2 Operation schematic for establishing the artery-vein track

12.9 Choose the suitable delivery system (mainly consisting of delivery sheath, dilator, loader, pusher and etc.) for the occluder, following the guidelines provided in Table 1.

Prepare the delivery system according to the instructions outlined in its IFU.

12.10 Insert the guide wire into the delivery sheath and guide it into the left ventricle.

Select an appropriate-sized occluder based on the results of angiography; typically, the occluder's disc diameter should be 2 mm larger than the determined diameter.

12.11 Place the occluder in heparinized saline. Pass the pusher through the loader and screw the device to the tip of the pusher by rotating it clockwise for four turns. After it is firmly attached, immerse the device and loader in the sterile heparinized normal saline (HepNS) and slowly pull the device into the loader. Use a syringe to inject heparin saline through the side arm of the delivery sheath to remove any bubbles.

12.12 Flush the delivery sheath and dilator with heparin saline. Insert the dilator into the delivery sheath and secure it with the locking mechanism. Guide the dilator/delivery sheath assembly over the guidewire into the left ventricle, following echo or X-ray guidance. Carefully withdraw the dilator and guidewire, leaving the delivery sheath in place within the left atrium.

12.13 Insert the loader into the longer sheath and advance it under X-ray guidance.

Open the disc near the left ventricle and pull it back toward the ventricular septum. If resistance is encountered around the ventricular septum as confirmed by echocardiography, secure the pusher, withdraw the sheath, and release the occluder disc near the right ventricle. Gently push and pull the occluder to ensure it is securely fixed. Verify with echocardiography that it does not interfere with the functioning of the aortic valve and tricuspid valve. Additionally, confirm through angiography that it does not affect the aortic valve's function. If both criteria are met, rotate the pushing pole counterclockwise to release the occluder.




Fig. 3 Operation schematic for VSD

12.14 During the operation, TTE is conducted using the five-chamber view at the cardiac apex and the short-axis view. The occluder should only be released when both views confirm that it is positioned closely adjacent to the residual edge of the ventricular septum.

12.15 Slowly retract the distal end of the pusher into the delivery sheath and withdraw the sheath. Complete the operation.

12.16 After the operation, the patient should be administered 3–5 mg/kg aspirin orally for 6 months.

### 13. Disposal



Treat all disposable devices appropriately according to the local requirements for medical device waste disposal.

### 14. Storage and Expiry Date












14.1 Store in a dry place, in a well-ventilated room with a relative humidity of no more than 80%, and with no corrosive gases or ultraviolet radiation.

14.2 This product has an expiry date of five years from the date of manufacture of this product if it stored under specified condition.

### 15. Symbol Definitions

SYMBOL	EXPLANATION OF SYMBOL
	MANUFACTURER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY



	DATE OF MANUFACTURE
	USE-BY DATE
<b>LOT</b>	BATCH CODE
<b>REF</b>	CATALOGUE NUMBER
<b>SN</b>	SERIAL NUMBER
<b>STERILE EO</b>	STERILIZED USING ETHYLENE OXIDE
	DOUBLE STERILE BARRIER SYSTEM AND STERILIZED USING ETHYLENE OXIDE
	DOUBLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE
	DO NOT RESTERILIZE
	DO NOT RE-USE
	DO NOT USE IF PACKAGE IS DAMAGED
	REFER TO INSTRUCTIONS FOR USE
<b>MD</b>	MEDICAL DEVICE
<b>UDI</b>	UNIQUE DEVICE IDENTIFIER
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	MR CONDITIONAL
<b>CE 2797</b>	CE MARKING AND IDENTIFICATION NUMBER OF NOTIFIED BODY
The color of the EO Sterilization indicator on the package bag turns yellow after EO sterilization.	

Shanghai Shape Memory Alloy Co., Ltd is not responsible for any direct, incidental, or consequential damages resulting from the abnormal use of the product.

### 17. Specification and Recommend sheath size

Table1 Device Specifications and Recommended Delivery System

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
SQFDQ- I a04	04	8.0±1.0	5.0±1.5	4.0±1.0	8.0±1.0	7-8F	Muscular defect
SQFDQ- I a05	05	9.0±1.0	5.0±1.5	5.0±1.0	9.0±1.0	7-8F	
SQFDQ- I a06	06	10.0±1.0	5.0±1.5	6.0±1.0	10.0±1.0	8-9F	
SQFDQ- I a07	07	11.0±1.0	5.0±1.5	7.0±1.0	11.0±1.0	8-9F	
SQFDQ- I a08	08	12.0±1.0	5.0±1.5	8.0±1.0	12.0±1.0	8-9F	
SQFDQ- I a09	09	13.0±1.0	5.0±1.8	9.0±1.0	13.0±1.0	8-9F	
SQFDQ- I a10	10	14.0±1.0	5.0±1.8	10.0±1.2	14.0±1.0	9-10F	
SQFDQ- I a12	12	16.0±1.0	5.0±1.8	12.0±1.2	16.0±1.0	9-10F	
SQFDQ- I a14	14	18.0±1.0	5.0±1.8	14.0±1.5	18.0±1.0	10-12F	
SQFDQ- I a16	16	20.0±1.0	5.0±1.8	16.0±1.5	20.0±1.0	10-12F	
SQFDQ- I a18	18	22.0±1.0	5.0±1.8	18.0±1.5	22.0±1.0	10-12F	
SQFDQ- I b04	04	10.0±1.0	7.0±1.5	4.0±1.0	8.0±1.0	7-8F	
SQFDQ- I b05	05	11.0±1.0	7.0±1.5	5.0±1.0	9.0±1.0	7-8F	
SQFDQ- I b06	06	12.0±1.0	7.0±1.5	6.0±1.0	10.0±1.0	8-9F	
SQFDQ- I b07	07	13.0±1.0	7.0±1.5	7.0±1.0	11.0±1.0	8-9F	
SQFDQ- I b08	08	14.0±1.0	7.0±1.5	8.0±1.0	12.0±1.0	8-9F	Muscular defect
SQFDQ- I b09	09	15.0±1.0	7.0±1.8	9.0±1.0	13.0±1.0	8-9F	
SQFDQ- I b10	10	16.0±1.0	7.0±1.8	10.0±1.2	14.0±1.0	9-10F	

### 16. DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

Descriptions or specifications in any content provided by Shanghai Shape Memory Alloy Co., Ltd., including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
SQFDQ- I b12	12	18.0±1.0	7.0±1.8	12.0±1.2	16.0±1.0	9-10F	
SQFDQ- I b14	14	20.0±1.0	7.0±1.8	14.0±1.5	18.0±1.0	10-12F	
SQFDQ- I b16	16	22.0±1.0	7.0±1.8	16.0±1.5	20.0±1.0	10-12F	
SQFDQ- I b18	18	24.0±1.0	7.0±1.8	18.0±1.5	22.0±1.0	10-12F	
SQFDQ- I c04	04	14.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	
SQFDQ- I c05	05	15.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
SQFDQ- I c06	06	16.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
SQFDQ- I c07	07	17.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
SQFDQ- I c08	08	18.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
SQFDQ- I c09	09	19.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	
SQFDQ- I c10	10	20.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
SQFDQ- I c12	12	22.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
SQFDQ- I c14	14	24.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	
SQFDQ- I c16	16	26.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	
SQFDQ- I c18	18	28.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	
SQFDQ- I d04	04	18.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	Muscular defect
SQFDQ- I d05	05	19.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
SQFDQ- I d06	06	20.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
SQFDQ- I d07	07	21.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
SQFDQ- I d08	08	22.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
SQFDQ- I d09	09	23.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
SQFDQ- I d10	10	24.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
SQFDQ- I d12	12	26.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
SQFDQ- I d14	14	28.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	
SQFDQ- I d16	16	30.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	
SQFDQ- I d18	18	32.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	
WTSQFDQ- I a04	04	8.0±1.0	5.0±1.5	4.0±1.0	8.0±1.0	7-8F	Muscular defect
WTSQFDQ- I a05	05	9.0±1.0	5.0±1.5	5.0±1.0	9.0±1.0	7-8F	
WTSQFDQ- I a06	06	10.0±1.0	5.0±1.5	6.0±1.0	10.0±1.0	8-9F	
WTSQFDQ- I a07	07	11.0±1.0	5.0±1.5	7.0±1.0	11.0±1.0	8-9F	
WTSQFDQ- I a08	08	12.0±1.0	5.0±1.5	8.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I a09	09	13.0±1.0	5.0±1.8	9.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I a10	10	14.0±1.0	5.0±1.8	10.0±1.2	14.0±1.0	9-10F	
WTSQFDQ- I a12	12	16.0±1.0	5.0±1.8	12.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I a14	14	18.0±1.0	5.0±1.8	14.0±1.5	18.0±1.0	10-12F	
WTSQFDQ- I a16	16	20.0±1.0	5.0±1.8	16.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I a18	18	22.0±1.0	5.0±1.8	18.0±1.5	22.0±1.0	10-12F	
WTSQFDQ- I b04	04	10.0±1.0	7.0±1.5	4.0±1.0	8.0±1.0	7-8F	
WTSQFDQ- I b05	05	11.0±1.0	7.0±1.5	5.0±1.0	9.0±1.0	7-8F	
WTSQFDQ- I b06	06	12.0±1.0	7.0±1.5	6.0±1.0	10.0±1.0	8-9F	
WTSQFDQ- I b07	07	13.0±1.0	7.0±1.5	7.0±1.0	11.0±1.0	8-9F	
WTSQFDQ- I b08	08	14.0±1.0	7.0±1.5	8.0±1.0	12.0±1.0	8-9F	

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
WTSQFDQ- I b09	09	15.0±1.0	7.0±1.8	9.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I b10	10	16.0±1.0	7.0±1.8	10.0±1.2	14.0±1.0	9-10F	
WTSQFDQ- I b12	12	18.0±1.0	7.0±1.8	12.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I b14	14	20.0±1.0	7.0±1.8	14.0±1.5	18.0±1.0	10-12F	
WTSQFDQ- I b16	16	22.0±1.0	7.0±1.8	16.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I b18	18	24.0±1.0	7.0±1.8	18.0±1.5	22.0±1.0	10-12F	
WTSQFDQ- I c04	04	14.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	
WTSQFDQ- I c05	05	15.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
WTSQFDQ- I c06	06	16.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I c07	07	17.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	
WTSQFDQ- I c08	08	18.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	
WTSQFDQ- I c09	09	19.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F	
WTSQFDQ- I c10	10	20.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F	
WTSQFDQ- I c12	12	22.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F	
WTSQFDQ- I c14	14	24.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F	
WTSQFDQ- I c16	16	26.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F	
WTSQFDQ- I c18	18	28.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F	
WTSQFDQ- I d04	04	18.0±1.0	10.0±1.5	4.0±1.0	10.0±1.0	7-8F	Muscular defect
WTSQFDQ- I d05	05	19.0±1.0	10.0±1.5	5.0±1.0	11.0±1.0	7-8F	
WTSQFDQ- I d06	06	20.0±1.0	10.0±1.5	6.0±1.0	12.0±1.0	8-9F	
WTSQFDQ- I d07	07	21.0±1.0	10.0±1.5	7.0±1.0	13.0±1.0	8-9F	

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect	
WTSQFDQ- I d08	08	22.0±1.0	10.0±1.5	8.0±1.0	14.0±1.0	8-9F	Muscular defect	
WTSQFDQ- I d09	09	23.0±1.0	10.0±1.8	9.0±1.0	15.0±1.0	8-9F		
WTSQFDQ- I d10	10	24.0±1.0	10.0±1.8	10.0±1.2	16.0±1.0	9-10F		
WTSQFDQ- I d12	12	26.0±1.0	10.0±1.8	12.0±1.2	18.0±1.0	9-10F		
WTSQFDQ- I d14	14	28.0±1.0	10.0±1.8	14.0±1.5	20.0±1.0	10-12F		
WTSQFDQ- I d16	16	30.0±1.0	10.0±1.8	16.0±1.5	22.0±1.0	10-12F		
WTSQFDQ- I d18	18	32.0±1.0	10.0±1.8	18.0±1.5	24.0±1.0	10-12F		
SQFDQ-IIa04	04	8.0±0.8	1.8±0.5	4.0±0.8	8.0±0.8	6-7F		Membranous defect
SQFDQ-IIa05	05	9.0±0.8	1.8±0.5	5.0±0.8	9.0±0.8	6-7F		
SQFDQ-IIa06	06	10.0±0.8	1.8±0.5	6.0±0.8	10.0±0.8	7-8F		
SQFDQ-IIa07	07	11.0±0.8	1.8±0.5	7.0±0.8	11.0±0.8	7-8F		
SQFDQ-IIa08	08	12.0±0.8	1.8±0.5	8.0±0.8	12.0±0.8	7-8F		
SQFDQ-IIa09	09	13.0±0.8	1.8±0.5	9.0±0.8	13.0±0.8	8-9F		
SQFDQ-IIa10	10	14.0±0.8	1.8±0.5	10.0±0.8	14.0±0.8	8-9F		
SQFDQ-IIa12	12	16.0±0.8	1.8±0.5	12.0±0.8	16.0±0.8	9-10F		
SQFDQ-IIa14	14	18.0±0.8	1.8±0.5	14.0±0.8	18.0±0.8	9-10F		
SQFDQ-IIa16	16	20.0±0.8	1.8±0.5	16.0±0.8	20.0±0.8	10-12F		
SQFDQ-IIa18	18	24.0±0.8	1.8±0.5	18.0±0.8	22.0±0.8	10-12F		
SQFDQ-IIa20	20	26.0±0.8	1.8±0.5	20.0±0.8	24.0±0.8	12-14F		
SQFDQ-IIb04	04	8.0±0.8	3.5±1.0	4.0±0.8	8.0±0.8	6-7F		
SQFDQ-IIb05	05	9.0±0.8	4.0±1.0	5.0±0.8	9.0±0.8	6-7F		



Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
SQFDQ-IIb06	06	10.0±0.8	4.0±1.0	6.0±1.0	10.0±0.8	7-8F	Membranous defect
SQFDQ-IIb07	07	11.0±1.0	4.0±1.0	7.0±1.0	11.0±1.0	7-8F	
SQFDQ-IIb08	08	12.0±1.0	4.0±1.0	8.0±1.0	12.0±1.0	7-8F	
SQFDQ-IIb09	09	13.0±1.0	4.5±1.0	9.0±1.2	13.0±1.0	8-9F	
SQFDQ-IIb10	10	14.0±1.5	4.5±1.0	10.0±1.2	14.0±1.5	8-9F	
SQFDQ-IIb12	12	16.0±1.5	4.5±1.0	12.0±1.5	15.0±1.5	9-10F	
SQFDQ-IIb14	14	18.0±1.5	4.5±1.0	14.0±1.5	17.0±1.5	9-10F	
SQFDQ-IIb16	16	22.0±1.5	5.0±1.0	16.0±1.5	20.0±1.5	10-12F	
SQFDQ-IIb18	18	24.0±1.5	5.0±1.0	18.0±1.8	22.0±1.5	10-12F	
SQFDQ-IIb20	20	26.0±1.5	5.0±1.0	20.0±1.8	24.0±1.5	12-14F	
WTSQFDQ-IIa04	04	8.0±0.8	1.8±0.5	4.0±0.8	8.0±0.8	6-7F	
WTSQFDQ-IIa05	05	9.0±0.8	1.8±0.5	5.0±0.8	9.0±0.8	6-7F	
WTSQFDQ-IIa06	06	10.0±0.8	1.8±0.5	6.0±0.8	10.0±0.8	7-8F	
WTSQFDQ-IIa07	07	11.0±0.8	1.8±0.5	7.0±0.8	11.0±0.8	7-8F	
WTSQFDQ-IIa08	08	12.0±0.8	1.8±0.5	8.0±0.8	12.0±0.8	7-8F	
WTSQFDQ-IIa09	09	13.0±0.8	1.8±0.5	9.0±0.8	13.0±0.8	8-9F	
WTSQFDQ-IIa10	10	14.0±0.8	1.8±0.5	10.0±0.8	14.0±0.8	8-9F	
WTSQFDQ-IIa12	12	16.0±0.8	1.8±0.5	12.0±0.8	16.0±0.8	9-10F	
WTSQFDQ-IIa14	14	18.0±0.8	1.8±0.5	14.0±0.8	18.0±0.8	9-10F	
WTSQFDQ-IIa16	16	20.0±0.8	1.8±0.5	16.0±0.8	20.0±0.8	10-12F	
WTSQFDQ-IIa18	18	24.0±0.8	1.8±0.5	18.0±0.8	22.0±0.8	10-12F	

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
WTSQFDQ-IIa20	20	26.0±0.8	1.8±0.5	20.0±0.8	24.0±0.8	12-14F	Membranous defect
WTSQFDQ-IIb04	04	8.0±0.8	3.5±1.0	4.0±0.8	8.0±0.8	6-7F	
WTSQFDQ-IIb05	05	9.0±0.8	4.0±1.0	5.0±0.8	9.0±0.8	6-7F	
WTSQFDQ-IIb06	06	10.0±0.8	4.0±1.0	6.0±1.0	10.0±0.8	7-8F	
WTSQFDQ-IIb07	07	11.0±1.0	4.0±1.0	7.0±1.0	11.0±1.0	7-8F	
WTSQFDQ-IIb08	08	12.0±1.0	4.0±1.0	8.0±1.0	12.0±1.0	7-8F	
WTSQFDQ-IIb09	09	13.0±1.0	4.5±1.0	9.0±1.2	13.0±1.0	8-9F	
WTSQFDQ-IIb10	10	14.0±1.5	4.5±1.0	10.0±1.2	14.0±1.5	8-9F	
WTSQFDQ-IIb12	12	16.0±1.5	4.5±1.0	12.0±1.5	15.0±1.5	9-10F	
WTSQFDQ-IIb14	14	18.0±1.5	4.5±1.0	14.0±1.5	17.0±1.5	9-10F	
WTSQFDQ-IIb16	16	22.0±1.5	5.0±1.0	16.0±1.5	20.0±1.5	10-12F	
WTSQFDQ-IIb18	18	24.0±1.5	5.0±1.0	18.0±1.8	22.0±1.5	10-12F	
WTSQFDQ-IIb20	20	26.0±1.5	5.0±1.0	20.0±1.8	24.0±1.5	12-14F	
SQFDQ-III04	04	12.0±1.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
SQFDQ-III05	05	13.0±1.0	4.0±1.5	5.0±1.0	9.0±1.0	8-9F	
SQFDQ-III06	06	14.0±1.0	4.0±1.5	6.0±1.0	10.0±1.0	8-9F	
SQFDQ-III07	07	15.0±1.0	4.0±1.5	7.0±1.2	11.0±1.0	8-9F	
SQFDQ-III08	08	16.0±1.2	4.0±1.5	8.0±1.2	12.0±1.2	9-10F	
SQFDQ-III09	09	17.0±1.2	4.5±1.5	9.0±1.2	13.0±1.2	9-10F	
SQFDQ-III10	10	18.0±1.2	4.5±1.5	10.0±1.5	14.0±1.2	9-10F	
SQFDQ-III12	12	20.0±1.5	4.5±1.5	12.0±1.5	16.0±1.2	10-12F	

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
SQFDQ-III14	14	22.0±1.5	4.5±1.5	14.0±1.8	18.0±1.5	10-12F	
SQFDQ-III16	16	24.0±1.5	5.0±1.5	16.0±1.8	20.0±1.5	10-12F	
SQFDQ-III18	18	26.0±1.5	5.0±1.5	18.0±1.8	22.0±1.5	12-14F	
WTSQFDQ-III04	04	12.0±1.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
WTSQFDQ-III05	05	13.0±1.0	4.0±1.5	5.0±1.0	9.0±1.0	8-9F	
WTSQFDQ-III06	06	14.0±1.0	4.0±1.5	6.0±1.0	10.0±1.0	8-9F	
WTSQFDQ-III07	07	15.0±1.0	4.0±1.5	7.0±1.2	11.0±1.0	8-9F	
WTSQFDQ-III08	08	16.0±1.2	4.0±1.5	8.0±1.2	12.0±1.2	9-10F	
WTSQFDQ-III09	09	17.0±1.2	4.5±1.5	9.0±1.2	13.0±1.2	9-10F	
WTSQFDQ-III10	10	18.0±1.2	4.5±1.5	10.0±1.5	14.0±1.2	9-10F	
WTSQFDQ-III12	12	20.0±1.5	4.5±1.5	12.0±1.5	16.0±1.2	10-12F	Membranous defect
WTSQFDQ-III14	14	22.0±1.5	4.5±1.5	14.0±1.8	18.0±1.5	10-12F	
WTSQFDQ-III16	16	24.0±1.5	5.0±1.5	16.0±1.8	20.0±1.5	10-12F	
WTSQFDQ-III18	18	26.0±1.5	5.0±1.5	18.0±1.8	22.0±1.5	12-14F	
SQFDQ-IV04	04	9.0±2.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
SQFDQ-IV05	05	10.0±2.0	3.5±1.5	5.0±0.8	9.0±1.0	7-8F	
SQFDQ-IV06	06	11.0±2.0	4.0±1.5	6.0±1.0	10.0±1.0	7-8F	
SQFDQ-IV07	07	12.0±2.5	4.0±1.5	7.0±1.0	11.0±1.0	8-9F	
SQFDQ-IV08	08	13.0±2.5	4.5±1.5	8.0±1.2	12.0±1.2	8-9F	
SQFDQ-IV09	09	14.0±2.5	5.0±1.5	9.0±1.2	13.0±1.2	9-10F	Membranous defect
SQFDQ-IV10	10	17.0±2.5	5.0±1.5	10.0±1.7	15.0±1.2	9-10F	

Catalog No.	Device Size	A LV Disc Diameter /mm	H Height of connecting waist /mm	B Connecting Waist diameter /mm	C RV Disc Diameter /mm	Smallest Recommended Delivery System Size	Applicable defect
SQFDQ-IV12	12	20.0±3.5	5.0±1.5	12.0±1.7	18.0±1.2	10-12F	
SQFDQ-IV14	14	22.0±3.5	5.0±1.5	14.0±1.8	20.0±1.5	10-12F	
SQFDQ-IV16	16	24.0±3.5	5.0±1.5	16.0±1.8	22.0±1.5	10-12F	
WTSQFDQ-IV04	04	9.0±2.0	3.5±1.5	4.0±0.8	8.0±1.0	7-8F	
WTSQFDQ-IV05	05	10.0±2.0	3.5±1.5	5.0±0.8	9.0±1.0	7-8F	
WTSQFDQ-IV06	06	11.0±2.0	4.0±1.5	6.0±1.0	10.0±1.0	7-8F	
WTSQFDQ-IV07	07	12.0±2.5	4.0±1.5	7.0±1.0	11.0±1.0	8-9F	
WTSQFDQ-IV08	08	13.0±2.5	4.5±1.5	8.0±1.2	12.0±1.2	8-9F	
WTSQFDQ-IV09	09	14.0±2.5	5.0±1.5	9.0±1.2	13.0±1.2	9-10F	
WTSQFDQ-IV10	10	17.0±2.5	5.0±1.5	10.0±1.7	15.0±1.2	9-10F	
WTSQFDQ-IV12	12	20.0±3.5	5.0±1.5	12.0±1.7	18.0±1.2	10-12F	
WTSQFDQ-IV14	14	22.0±3.5	5.0±1.5	14.0±1.8	20.0±1.5	10-12F	
WTSQFDQ-IV16	16	24.0±3.5	5.0±1.5	16.0±1.8	22.0±1.5	10-12F	Membranous defect

Note: Device applicable for muscular ventricular septal defect cannot be used for membranous ventricular septal defect and vice versa.

**18.E-IFU**

The exact same e-IFU (pdf version) can also be found on the website: [https:// en.scientechmed.com/product/](https://en.scientechmed.com/product/).

Note:Whenever the manufacturer's instructions for use are updated, the revised version will be promptly uploaded. Due to the challenge of individually notifying every end user about these changes, we recommend that customers regularly browse and check the website for updates.

In addition, a computer installed windows system, one currently common browser and pdf file reader is required.